


b e n t e c

**STERILE**



## Silicone Elastomer Tubing

- Multiple tubing sizes from 2 to 28 French.
- Smooth, soft, and pliable.
- Tubing is available in translucent and radiopaque.
- Custom lengths by request.



To order, call 800•767•9175  
Visa/Mastercard/Procurement Cards Accepted

1380 E. Beamer St.  
Woodland, CA 95776

Fax: 888-828-8869  
Internet: [www.bentecmed.com](http://www.bentecmed.com)

## Bentec Silicone Elastomer Tubing

STERILE

Manufactured from the highest quality Silastic™ Silicone, Bentec Silicone Elastomer Tubing is smooth, soft and pliable, making it easier to use and more comfortable for the patient. Provided **sterile** in translucent and radiopaque. Available in a wide variety of French sizes, in 100 cm. (39.37 in.) length. Latex-free for a safer standard of care.

100 cm. length (39.37 in.)

Size	Dimensions		Translucent	Radiopaque	Size	Dimensions		Translucent	Radiopaque
2 Fr.	OD (.025)	ID (.012)	GS75160-02	GS75160-42	14 Fr.	OD (.183)	ID (.122)	GS75160-14	GS75160-54
3 Fr.	OD (.037)	ID (.020)	GS75160-03	GS75160-43	15 Fr.	OD (.197)	ID (.113)	GS75160-15	GS75160-55
4 Fr.	OD (.052)	ID (.027)	GS75160-04	GS75160-44	16 Fr.	OD (.210)	ID (.140)	GS75160-16	GS75160-56
5 Fr.	OD (.065)	ID (.030)	GS75160-05	GS75160-45	18 Fr.	OD (.240)	ID (.160)	GS75160-18	GS75160-58
6 Fr.	OD (.078)	ID (.047)	GS75160-06	GS75160-46	19 Fr.	OD (.250)	ID (.113)	GS75160-19	GS75160-59
7 Fr.	OD (.092)	ID (.055)	GS75160-07	GS75160-47	20 Fr.	OD (.262)	ID (.174)	GS75160-20	GS75160-60
8 Fr.	OD (.105)	ID (.063)	GS75160-08	GS75160-48	22 Fr.	OD (.290)	ID (.200)	GS75160-22	GS75160-62
9 Fr.	OD (.118)	ID (.071)	GS75160-09	GS75160-49	24 Fr.	OD (.314)	ID (.209)	GS75160-24	GS75160-64
10 Fr.	OD (.131)	ID (.078)	GS75160-10	GS75160-50	26 Fr.	OD (.335)	ID (.295)	GS75160-26	
11 Fr.	OD (.142)	ID (.110)	GS75160-11	GS75160-51	28 Fr.	OD (.360)	ID (.181)	GS75160-28	
12 Fr.	OD (.155)	ID (.105)	GS75160-12	GS75160-52					

ID = Inner Diameter OD = Outer Diameter All dimensions provided in inches

## Bentec Medical Grade Silicone Sheeting

STERILE

Available either mesh reinforced or non-reinforced for a wide variety of applications. Provided **sterile**. Additional sheeting sizes up to 14" x 18" available.

Non-Reinforced	Thickness	Size	Qty. per pkg.	Reinforced	Thickness	Size	Qty. per pkg.
PR72034-20N	.005	3 x 4	2	PR72034-007R	.007	6 x 8	1
PR72034-00N	.005	6 x 8	1	PR72034-01R	.010	6 x 8	1
PR72034-21N	.010	3 x 4	2	PR72034-02R	.020	6 x 8	1
PR72034-01N	.010	6 x 8	1	PR72034-03R	.030	6 x 8	1
PR72034-22N	.020	3 x 4	2	PR72034-04R	.040	6 x 8	1
PR72034-02N	.020	6 x 8	1	PR72034-05R	.050	6 x 8	1
PR72034-23N	.030	3 x 4	2	PR72034-06R	.060	6 x 8	1
PR72034-03N	.030	6 x 8	1	PR72034-07R	.070	6 x 8	1
PR72034-24N	.040	3 x 4	2	PR72034-08R	.080	6 x 8	1
PR72034-04N	.040	6 x 8	1	PR72034-10R	.100	6 x 8	1
PR72034-05N	.050	6 x 8	1	PR72034-11N	.120	6 x 8	1
PR72034-06N	.060	6 x 8	1	PR72034-12N	.125	6 x 8	1
PR72034-07N	.070	6 x 8	1				
PR72034-08N	.080	6 x 8	1				
PR72034-09N	.090	6 x 8	1				

All Bentec products are latex-free

## Bentec Silicone Block

STERILE

Block design allows the surgeon to custom fabricate, at surgery, an implantable prosthesis for a specific surgical indication. Provided **sterile** in four consistencies: soft, medium, firm and extra firm.

Block	Size	Consistency
PR88032-01	7cm x 12cm x 1.5cm	Soft
PR88032-02	"	Medium
PR88032-03	"	Firm
PR88032-04	"	X-Firm
PR88032-21	7cm x 4cm x 1.5cm	Soft
PR88032-22	"	Medium
PR88032-23	"	Firm
PR88032-24	"	X-Firm



Visa/Mastercard/Procurement  
Cards Accepted  
Special Orders Available

Internet: [www.bentecmed.com](http://www.bentecmed.com)  
E-mail: [bentec@bentecmed.com](mailto:bentec@bentecmed.com)

To order, call 800•767•9175